



CORNER SANIGEL CSD-600MLW

AUTOMATIC SANIGEL TOILET SEAT WIPES

An effective and quick-drying germicidal gel system used to clean toilet seats safely. Positioned inside toilet cubicles, the SANIGEL uses the same tissue paper normally provided in the toilets.

SANIGEL system is simple and hassle-free.

For convenience and hygiene purposes the SANIGEL DISPENSER has a sensor activator, so that it can be a touch-free operation. A push-bar activator is also built-in, if required.



AUTOMATIC SANIGEL ASDT-950S

USING A SANIGEL INVOLVES 4 SIMPLE STEPS:

1. Pull out approximately 2 feet length of toilet paper and fold it a few times.
2. Position the folded toilet paper under the dispenser and the SANIGEL will be automatically dispensed.
3. Clean the toilet seat with a wiping action.
4. Dispose off the used toilet paper by throwing into the toilet bowl.

SANIGEL IS A MILD AND EFFECTIVE ALCOHOL-BASE CLEANSER THAT IS SAFE FOR SKIN CONTACT OR HANDS..

MODELS/DIMENSIONS:

AUTOMATIC SANIGEL ASDT-950S: 24CM(H) X 14CM(W) X 9CM(D)
 MANUAL SANIGEL SDW-950S: 24CM(H) X 14CM(W) X 9CM(D)
 CORNER SANIGEL CSD-600MLW: 21CM(H) X 15.5CM(W) X 8.5CM(D)

REFILLS PACKING:

12 X 600ml disposable sachet per carton
 12 X 1000ml disposable sachet per carton

